

Acid Catalyzed Transfer Hydrogenation


Optimizing the Kinetics

Continuous-flow Catalytic Asymmetric Hydrogenations: Reaction Optimization Using FTIR Inline Analysis

The majority of enantioselective reductions depend on various transition-metal-catalyzed high-pressure hydrogenations. However, the drawbacks of these method conditions are often tolerated due to the excellent enantio selectivities and activity they obtain.


Rueping et. al studied the asymmetric organocatalytic hydrogenation of benzoxazines in continuous-flow microreactors as a metal-free alternative. Coupling inline mid-infrared spectroscopy for real-time reaction monitoring enabled the group to quickly, optimize the reaction for duration (kinetics) and yield (product quality) while minimizing the volume of material needed.

By studying the effects of a select set of parameters (temperature, flow rate, and residence time) the optimum set of conditions were determined as based on overall product yield (Figure 1). Inline mid-infrared provided the instantaneous feedback as each of these parameters were examined for their beneficial or negative contribution to the chemistry.



Entry	1a [mol %]	Residence Time [min]	Flow Rate [mL min ⁻¹]	Yield [%] ^P
1	2	20	0.1	50%
2	2	40	0.1	87%
3	2	60	0.1	98%

Figure 1. Optimization of the Bronsted acid catalyzed reduction of benzoxazines. Reaction Conditions: **3a**, **2a** (1.2 equiv), **1a** in CHCl_3 (0.05M) at 60 °C. ^PIsolated yields after column chromatography.








Figure 2. *In situ* ReactIR monitoring: (a) Trend curve of product formation at different temperatures. (b) Reaction spectra showing the consumption of the substrate and formation of the product at different temperatures. (c) Three-dimensional time-resolved spectral data.



Results

In situ conversion, as measured by ReactIR, tracks changes to the substrate and product over the course of the experiment. This allows researchers to determine key parameters for control and optimization. For example, by varying the temperature from 5 °C to 60 °C (5 °C increments), the ideal reaction temperature was identified to be 60 °C in a single eight-hour reaction. Additionally, ReactIR trends provided time-resolved reaction conversion data measuring the relative concentration of both product and substrate.

Conclusions


In this study, ReactIR demonstrated the ability to detect and monitor the substrate and product kinetic behavior leading to the optimized reactions conditions: temperature 60 °C, flow rate 0.1 mL/min, and 98 % yield (confirmed by column chromatography).

ReactIR: Real-Time Reaction Mechanisms




Monitor Chemistry *In Situ*

Monitor chemistry as it takes place (in flow) to achieve optimal reaction mechanics for a greener, safer process in the shortest time frame possible.



Comprehensive Reaction Information in Real Time

Detect and monitor instantaneous changes (initiation, steady state, process deviation, and endpoint) that are otherwise difficult or impossible to detect via offline methods.



Reaction Kinetics, Pathway and Mechanism

In-depth understanding of reaction mechanism for process efficiency, robustness and safety.



Applications of ReactIR include:

- Reaction Kinetics
- Flow Chemistry
- High Pressure Reactions
- Polymer Synthesis
- Energetic Reactions

Mettler-Toledo AutoChem, Inc.

7075 Samuel Morse Drive
Columbia, MD 21046 USA
Telephone +1 410 910 8500
Fax +1 410 910 8600

Email autochem@mt.com
Internet www.mt.com/autochem

Subject to technical changes
© 01/2017 Mettler-Toledo AutoChem, Inc.

www.mt.com/ReactIR

For more information